

Florida Integrated Science Center Sirenia Project

2201 N.W. 40th Terrace Gainesville, Florida 32605-3574 e-mail: Robert Bonde@usgs.gov

Phone: (352) 264-3555

Fax: (352) 374-8080

SIRENIAN BLOOD COLLECTION PROTOCOL FOR GENETIC ANALYSIS

Ship To: U.S.G.S. - Sirenia Project

Manatee Genetics Program 2201 NW 40th Terrace

Gainesville, Florida USA 32605-3574

Email: rbonde@usgs.gov

Permits: Federal FWS Endangered Species Research Permit MA-791721

CITES Importation Permit US-808447/9

PRESERVING BLOOD WITH LYSIS BUFFER

Materials:

- 1) Lysis buffer: 100 mM Tris-HCl, pH 8; 100 mM EDTA, pH 8; 10 mM NaCl; 1.0% SDS (weight/volume)
- 2) Syringe or vacutainer apparatus without anticoagulant agents
- 3) Screw-cap tubes or other appropriated sealed containers

Methods:

- 1) Collect 1-2 ml whole blood in sterile syringe or vacutainer
- 2) Add blood immediately to a labeled tube containing lysis buffer. Add approximately a 1:10 ratio of blood to lysis buffer
- 3) Invert tube gently several times to mix
- 4) Samples can be stored at room temperature or in a refrigerator for up to 1 year. Avoid extended exposure to direct sunlight or high temperatures
- 5) Ship samples as soon as possible by air freight or express courier

Note: SDS is sodium dodecyl sulfate (=sodium lauryl sulfate). Lysis buffer is nontoxic and can be stored for an extended period at room temperature.